

A Clinical Evaluation of an Automated Software Program (Certegra® P3T® PA) for Patient Specific Contrast Injection During Chest CTA to Exclude Pulmonary Embolism

Christopher R Deible MD, PHD¹, Jacob Alexander MD¹, Iclal Ocak MD¹, Maryam Ghadimi Mahani MD¹, John Kalafut BS, MS², Janet E Durick MD¹, Carl R Fuhrman MD¹, Darlene Frasher RN, MSN¹, Karen M Pealer BA, CCRC¹, Michael P. Federle MD¹, Joan M Lacomis MD¹ Department of Radiology, University of Pittsburgh Medical Center, Pittsburgh, PA¹, Medrad, Inc., Pittsburgh, PA²

INTRODUCTION

The role of CT angiography (CTA) has been well established in the detection of pulmonary thromboembolism. According to the PIOPEd II study, CTA overall sensitivity and specificity is 83% and 96% respectively.^{1,2} Optimizing contrast injection and scanning parameters has become of increased importance with the faster multidetector scanners to achieve diagnostic quality images.²⁻⁶

The purpose of our study is to assess if a prototype automated software program for patient specific contrast injection, Certegra® P3T® PA (Pulmonary Angiography), Medrad, Inc., Pittsburgh, PA, is comparable to or offers advantages over our site specific standard protocol used for chest CTA to exclude pulmonary embolism (PE).

METHODS AND MATERIALS

62 emergency department patients referred for chest CTA to exclude PE underwent informed consent for this study and were randomized to Certegra® P3T® PA versus Standard (control) groups. All had 18 gauge IV access, received the contrast agent Ioversol (350mg/ml iodine; Mallinckrodt, St Louis, MO) and were scanned on one 64 slice CT scanner (VCT, GE Healthcare; Milwaukee, WI) by selected technologists monitored by selected investigators (JL, CD, JA). Scan parameters: 0.625 mm collimation; 0.24 pitch; 9.6 mm/sec table speed; 350 msec rotation time; 120 kV; 280 - 550 mAs. Recorded patient parameters: height, weight, age, sex and heart rate.

CTA Standard Group: Test bolus = 20cc contrast/ 50cc saline flush @ 4cc/sec to time for main pulmonary artery (MPA). Scan delay = time to contrast peak MPA + 9 seconds. Scan bolus = 80cc contrast/ 50cc saline flush @ 4cc/sec.

CTA Certegra® P3T® PA Group: As a safety measure, a default max which allowed injection rate of 6cc/sec was pre-selected. Height, weight, age, sex, heart rate and scan duration were entered into Certegra® P3T® PA which generated test bolus parameters. From the test bolus, time to peak and peak density (HU) in the MPA were entered into Certegra® P3T® PA which generated scan bolus and scan delay parameters.


Data Collection/Analysis: Two readers, blinded to injection method, jointly measured density (HU) of main (MPA, RPA, LPA) and segmental pulmonary arteries (bilateral upper and lower lobes), and SVC. Three other blinded readers qualitatively scored scans compared to an "adequate" example for image quality to assess for PE, noting limitations: poor contrast, motion, quantum mottle, SVC streak and artifact. The mean and standard deviation for each group was calculated separately. Statistical analysis was performed with the Student's t-test, and the Wilcoxon rank sum test.

Figure 2: Certegra® P3T® PA Scans


BMI 17.5
HR 96 bmp
CO ~4.25 L/min
Contrast volume: 78 cc



BMI 24.9
HR 100 bpm
CO ~6.2 L/min
Contrast volume: 84 cc



BMI: 20.5
HR: 89 bpm
CO ~7.12 L/min
Contrast volume: 87 cc



Estimated cardiac output was computed using standard look-up tables

A Clinical Evaluation of an Automated Software Program (Certegra® P3T® PA) for Patient Specific Contrast Injection During Chest CTA to Exclude Pulmonary Embolism

Table 1: Comparison of Certegra® P3T® PA Enhancement (t-test)





Figure 3: Certegra® P3T® PA versus Standard



Contrast enhancement of pulmonary arteries in Certegra® P3T® PA exams vs Standard exams ($p < 0.01$ to $p < 0.005$, Wilcoxon rank sum)

RESULTS

- ◆ Certegra® P3T® PA population: 20 women, 11 men
 - Age: mean 43.7 yrs (range: 20-76)
 - Weight: mean 180 lbs +/- 59, median 172 lbs
 - BMI: mean 29 +/- 9, median 26
- ◆ Standard population: 22 women, 9 men
 - Age: mean 44.4 yrs (range: 19- 89)
 - Weight: mean 181 lbs +/- 75, median 185 lbs
 - BMI: mean 30 +/- 10, median 29
- ◆ Higher average image quality score of Certegra® P3T® PA exams (mean 4.2 +/- 0.8) vs. Standard exams (mean 3.6 +/- 1.2) ($p < 0.05$; Wilcoxon rank sum)
- ◆ Higher percentage of exams ranked as diagnostic without limitation (positive or negative for PE) in the Certegra® P3T® PA exams (100%) vs. Standard exams (73%*) ($p < 0.05$ Wilcoxon rank sum)
 - * Contrast related problems (63%) were most frequent factors cited for scan limitation
- ◆ Average contrast dose (test + scan) for Certegra® P3T® PA scans was higher (114mls +/- 12mls: range 76-152) vs Standard scans (100 mls) (z-score 6.6, $p < .001$, 1 sample z-test)

CONCLUSIONS

In this small study, the prototype Certegra® P3T® PA automated software program for patient specific contrast injection offered improved, more consistent contrast enhancement of the target pulmonary arteries over a variety of patient specific parameters. There was a higher percentage of diagnostic quality exams albeit at a slightly higher contrast dose than our site's standard protocol used to exclude PE.

REFERENCES

- Stein PD, Fowler SE, Goodman LR, et al. Multidetector computed tomography for acute pulmonary embolism. *N Engl J Med.* 2006; 354(22): 2317-27.
- Stein PD, Hull RD. Multidetector computed tomography for the diagnosis of acute pulmonary embolism. *Curr Opin Pulm Med.* 2007; 13(5): 384-8
- Raptopoulos V, Boiselle PM. Multi-detector row spiral CT pulmonary angiography: comparison with single-detector row spiral CT. *Radiology.* 2001; 221(3): 606-13.
- Yankelevitz DF, Shaham D, et al. Optimization of contrast delivery for pulmonary CT angiography. *Clin Imaging.* 1998 Nov-Dec; 22(6): 398-403.
- Tillie-Leblond I, Mastora I, Radenne F, et al. Risk of pulmonary embolism after a negative spiral CT angiogram in patients with pulmonary disease: 1-year clinical follow-up study. *Radiology.* 2002; 223 (2): 461-7.
- MacKenzie JD, Nazario-Larrieu J, Cai T, et al. Reduced-dose CT: effect on reader evaluation in detection of pulmonary embolism. *Am J Roentgenol.* 2007; 189(6): 1371-9.
- Jones SE, Wittram C. The indeterminate CT pulmonary angiogram: imaging characteristics and patient clinical outcome. *Radiology* 2005; 237; 329-337.




Bayer HealthCare

Radiology
Bayer HealthCare Medical Care

Bayer HealthCare LLC
100 Bayer Boulevard
P.O. Box 915
Whippany, NJ 07981

www.ri.bayer.com

 Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780 USA

Customer Service/Orders
1-800-633-7231

Customer Service Fax
1-412-767-4120

Bayer reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Please contact your authorized Bayer HealthCare representative for the most current information.

Bayer, the Bayer Cross, Medrad®, Certegra® and P3T® are trademarks of the Bayer group of companies.

© 2014 Bayer

RAD-INF-14-06087 November 2014