



OUR VALUE ADDS UP




Centargo Disposables



Saline Tag

- + Required when delivering saline from the same source to multiple patients
- + Once a saline fluid source is punctured, the Saline Tag must be attached for the duration of use
- + Catalog number: CENT-SALINE-TAG



Replacement Spike

- + Replace the spike with every fluid container change*
- + Catalog number: CENT-RS



Patient Line

- + Snap-in single hand design
- + Auto-primers upon insertion
- + Two check valves, recessed connection to reduce contamination risk*
- + Single use, 250 cm
- + Catalog number: CENT-PL

Ready for the next patient in less than 20 seconds. Centargo auto primes the patient line upon connection.

*Discard disposables in case of suspected fluid path contamination.

MEDRAD® Centargo
CT Injection System



Day Set

- + Labeled for twenty four (24) hours of use or 25 bottles of contrast media, whichever comes first
- + Pre-assembled Day Set includes three reservoirs of 200 mL each (one saline, two contrast) and three spikes
- + Recessed connection to reduce contamination risk
- + Catalog number: CENT-DS



Centargo auto-fills the Day Set at the end of each patient exam.



Centargo daily set-up in less than 2 minutes.

Ready for the next patient in less than 20 seconds.

Disposables Ordering Information

Item	Catalog #	Quantity
Centargo Day Set	CENT-DS	Box of 4
Centargo Patient Line	CENT-PL	Box of 50
Centargo Replacement Spikes	CENT-RS	Box of 50
Centargo Saline Tag	CENT-SALINE-TAG	Box of 50

Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

Bayer, the Bayer Cross, MEDRAD, MEDRAD Centargo, and Centargo are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

© 2025 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.



Bayer HealthCare LLC
100 Bayer Boulevard
P.O. Box 915
Whippany, NJ 07981
U.S.A.
Phone: +1-412-767-2400
+1-800-633-7231
Fax: +1-412-767-4120



Manufacturer
Imaxeon Pty Ltd
Unit 1, 38 – 46 South Street
Rydalmere, NSW 2116
Australia
Phone: +1-61-2-8845-4999
Fax: +1-61-2-8845-4998
Imaxeon Pty Ltd is a Bayer company

[Learn more at centargo.com](https://www.centargo.com)