



# More than 8,000 injection systems in the US are connected. Yours can be, too.\*

Help improve uptime, streamline service and strengthen security.



Cortenic™ Connectivity



## WHY CONNECT YOUR BAYER INJECTOR?

### + Stay Secure While Minimizing Disruption

Receive remote cybersecurity updates and patching aligned with hospital IT requirements.

### + Spot Issues Before They Escalate

Enable remote diagnostics and expert technical support to identify issues early and efficiently resolve them, helping minimize downtime and service visits.

### + Get Support with Push for Service

Initiate support directly with Push for Service, streamlining access to Bayer Service and reducing reliance on traditional call flows.

### + Enable Software and Future-Ready Workflows

Establish a foundation for activating workflow software while supporting efficient implementation and ongoing technical support.

### + Integrate Across Your Systems

Connect to Modality Worklist, PACS, RIS and Speech Recognition systems through remote connectivity—supporting more seamless workflows today and future expansion over time.

## GETTING CONNECTED IS SIMPLE

### Three steps to connect:

1. Gain approval from your IT team
2. Ensure an active network jack in the control room
3. Enable outbound firewall communication on port 443



**Connect Today**

**Contact your Bayer Representative or visit [go.bayer.com/cortenicconnectivity](https://go.bayer.com/cortenicconnectivity)**

\*Connectivity is included as an entitlement of purchase for MEDRAD® Centargo and MEDRAD® Stellant FLEX CT Injection Systems and MEDRAD® MRXperion MR Injection Systems. Other legacy injectors are not eligible.

Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

Bayer, the Bayer Cross, MEDRAD, MEDRAD Centargo, Centargo, Stellant FLEX, MEDRAD Stellant FLEX, MEDRAD MRXperion, MRXperion and Cortenic are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

© 2026 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer. PP-COR-US-0110-1 July 2026