



**Stylish ergonomics designed for function**

**Configurations for  
Siemens Medical Systems**

**medrad<sup>®</sup> Mark 7 Arterion<sup>®</sup>**  
Injection System

# Siemens Standard Configurations for Medrad® Mark 7 Arterion® Injection System

Medrad® Mark 7 Arterion® Injection System offers multiple configurations for maximum mounting flexibility

## Siemens Catalog No: BART700PEDL *Medrad® Mark 7 Arterion® Pedestal System*

- ◆ Medrad® Mark 7 Arterion® Base System: Injector Head with 10' head cable, Display Control Unit, Power Supply
- ◆ Integrated Pedestal Unit with Articulating Arm and Fulcrum
- ◆ 10' Display Cable
- ◆ 12' Coiled Hand Switch
- ◆ 16.5' Siemens Axiom Artis ISI Cable (Including Zee)
- ◆ Syringe Heat Maintainer
- ◆ Printed Operation Manual Medrad® Mark 7 Arterion® Injection System
- ◆ Multilingual Operation and Training Manual on CD
- ◆ Service Manual (English Only)
- ◆ Standard 1 year Warranty Paperwork



## Siemens Catalog No: B59941325 *Medrad® Mark 7 Arterion® Integrated Pedestal for Use on the Siemens Cios Alpha C-arm System*

- ◆ Medrad® Mark 7 Arterion® Base System: Injector Head with 10' head cable, Display Control Unit, Power Supply\*
- ◆ Integrated Pedestal Unit with Articulating Arm and Fulcrum
- ◆ 10' Display Cable
- ◆ 12' Coiled Hand Switch
- ◆ 16.4' Siemens Cios Alpha Equipotential Cable
- ◆ Syringe Heat Maintainer
- ◆ Printed Operation Manual Medrad® Mark 7 Arterion® Injection System
- ◆ Multilingual Operation and Training Manual on CD
- ◆ Service Manual (English Only)
- ◆ Standard 1 year Warranty Paperwork



\*Installation of the Medrad® Mark 7 Arterion® Injection System is included with the Cios Alpha C-arm System

## Siemens Catalog No: BART700TABL *Medrad® Mark 7 Arterion® Rack Mount Package*

- ◆ Medrad® Mark 7 Arterion® Base System: Injector Head with 10' head cable, Display Control Unit, Power Supply
- ◆ Display Control Unit Desk Stand Kit
- ◆ 6' Coiled Hand Switch
- ◆ 90' Head Extension Cable
- ◆ 100' Display Control Unit Cable
- ◆ Cable Floor Bracket
- ◆ 16.5' Siemens Axiom Artis ISI Cable (Including Zee)
- ◆ Syringe Heat Maintainer
- ◆ Printed Operation Manual Medrad® Mark 7 Arterion® Injection System
- ◆ Multilingual Operation and Training Manual on CD
- ◆ Service Manual (English Only)
- ◆ Standard 1 year Warranty Paperwork



## Siemens Catalog No: BART700OCS *Medrad® Mark 7 Arterion® Overhead Counterpoise System: Injector Head*

- ◆ Injector head with 65' Head Extension Cable
- ◆ Display Control Unit, Power Supply
- ◆ Display Control Unit Desk Stand Kit
- ◆ 6' Coiled Hand Switch
- ◆ 75' Display Control Unit Cable
- ◆ 16.5' Siemens Axiom Artis ISI Cable (Including Zee)
- ◆ Syringe Heat Maintainer
- ◆ Printed Operations Manual Medrad® Mark 7 Arterion® Injection System
- ◆ Multilingual Operation and Training Manual on CD
- ◆ Service Manual (English Only)
- ◆ Standard 1 year Warranty Paperwork



## Rack and Pedestal Mount Package Mounting Options

**Siemens Catalog No: BKMA320RT**  
Free Standing Pedestal on Wheels

## Overhead Counterpoise System Package Mounting Options

**Siemens Catalog No: B3016427**  
OCS 580mm Trolley Mount, Portegra 2

**Siemens Catalog No: B3016425**  
OCS 800mm Trolley Mount, Portegra 2

**Siemens Catalog No: B3016424**  
OCS 580mm Ceiling Mount, Portegra 2

**Siemens Catalog No: B3016426**  
OCS 850mm Ceiling Mount, Portegra 2

## Accessories

**Siemens Catalog No: BART700VFL**  
Medrad® VFlow™ Software

**Siemens Catalog No: B59893002**  
Arterion Footswitch

**Siemens Catalog No: B59898551**  
Arterion Pressure Jacket

**Siemens Catalog No: BART700DCUDM**  
Desk Stand for Display

**Siemens Catalog No: BART700DCUTMA**  
Adjustable Table Mount for Display

**Siemens Catalog No: BART700DCUWM**  
Wall Mount for Display



\*As compared with Medrad® Mark V ProVis® Angiographic Injection System

Medrad® Mark 7 Arterion® has stylish ergonomics that are designed for function; the pedestal with a small footprint and reduced weight increases mobility around your busy lab. And the light head now has a handle to make it easier to position for injection.

## Variable Flow

Medrad® VFlow™ Sterile Hand Controller\* provides the flexibility to deliver hand injection with the accuracy of an injector to enable consistent contrast delivery.

<p><b>Control</b></p> <p>Real-time control of flow rate to enhance contrast delivery across procedure type</p>	<p><b>Consistency</b></p> <p>Automated injection for accurate and consistent contrast delivery every time</p>
<p><b>Flexibility</b></p> <p>Flexibility to deliver contrast in variable or fixed flow rates with a single injector system</p>	<p><b>Operator Safety</b></p> <p>Ability to control the injection procedure from a distance, reducing operator radiation exposure</p>



\* Refer to the Operators Manual for a description of the intended use, contraindications, warnings and instructions for use.

The Medrad® Mark 7 Arterion® Injection System is marketed for the administration of intravascular contrast media and common flushing agents during angiographic procedures. See Product Information for Use for specific and complete prescribing information.

# Equipment Service by Bayer

The #1 way to keep Bayer devices performing at peak efficiency is with Equipment Service by Bayer at the time of purchase.

**Think ahead** with VirtualCare®, our round-the-clock remote support.

**Plan ahead** with our proprietary predictive maintenance and service programs.

**Stay ahead** with our rapid response teams.



## Bayer

Pharmaceuticals Division

Bayer HealthCare LLC  
100 Bayer Boulevard  
P.O. Box 915  
Whippany, NJ 07981  
[www.radiologysolutions.bayer.com](http://www.radiologysolutions.bayer.com)



Bayer Medical Care Inc.  
1 Bayer Drive  
Indianola, PA 15051 USA  
Customer Service/Orders  
1-800-633-7231  
Customer Service Fax  
1-412-767-4120

Bayer reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

Bayer, the Bayer Cross, Medrad®, Mark 7 Arterion®, VFlow™, and Mark V ProVis® are trademarks of the Bayer group of companies. All other trademarks are properties of their respective owners.

© 2016 Bayer

PP-ART-US-0058 October 2016