



FLEX Solution Plans

POWERED BY CORTENIC™ CONNECTIVITY

Optimizing the capabilities and performance
of your connected MEDRAD® Stellant FLEX



MEDRAD Stellant FLEX

POWERED BY CORTENIC CONNECTIVITY

A Connected CT Solution for Various Suites and Configurations

MEDRAD Stellant FLEX, powered by Cortenic Connectivity, is a single-patient, dual syringe injection system designed to help you meet your financial, clinical and organizational goals.




Versatility, Enhanced.

An easy-to-use injector that adapts to your department's operational needs across a variety of CT suite workflows, so you can stay focused on patient care.



Indication, Expanded.

The first and only CT injection system currently 510(k) cleared for use with CEM, that enables consistent contrast delivery for enhanced visualization of breast lesions.¹



Performance, Maximized.

Technology built on proven MEDRAD Stellant reliability and backed by expert clinical and technical services.



Connectivity, Elevated.

Powered by Cortenic Connectivity, MEDRAD Stellant FLEX is a comprehensive solution to help drive performance across your connected suite, with integrated software, services, education and data.



Contrast, Optimized.

Streamline preparation and optimize in-suite operations with MEDRAD Stellant FLEX and the Imaging Bulk Package Transfer Set.

Solution Plans, Activated.

Optimize the capabilities and performance of your connected MEDRAD Stellant FLEX.

INTRODUCING CORTENIC CONNECTIVITY:

Cortenic Connectivity powers connected solutions from Bayer. Connectivity unlocks access to advanced features, which can do more to drive workflow efficiency, optimize performance and deliver data for smarter insights in the radiology suite.

Tiered Solution Plans

To unlock the advantages of connected software and services, pair MEDRAD Stellant FLEX with the solution plan that is right for your department. FLEX Solution Plans powered by Cortenic Connectivity combine Bayer's injector technology, software and services to activate integrated, scalable solutions that adapt to the needs of your CT suite.

FLEX Solution Plans elevate the features and benefits of the connected injector, and are designed to drive outcomes in your department, from results that streamline your essential daily operations, to those that drive comprehensive clinical and operational initiatives. Solution Plans bring the best of Bayer's technical and clinical knowledge to your suite, allowing the experts you trust to become an extension of your team.

Click to Explore
Our Plans

BASE TIER



ESSENTIAL SERVICE TIER



ESSENTIAL SOFTWARE TIER



COMPREHENSIVE TIER



A Tiered Approach to Connected Solutions



All customers begin with the Base Solution – a connected MEDRAD Stellant FLEX.

From there, the Essential Solution Plans offer the fundamentals in support and maintenance services or software.

The Comprehensive Solution Plan offers integrated software, advanced services, education and data for the connected suite.

BASE SOLUTION TIER

ESSENTIAL SERVICE SOLUTION TIER

ESSENTIAL SOFTWARE SOLUTION TIER

COMPREHENSIVE SOLUTION TIER

EXPLORE ALL TIER OPTIONS >>>

Get Connected with the MEDRAD Stellant FLEX Base Solution

Every customer starts with **The Base Solution** — a connected **MEDRAD Stellant FLEX**. When remote connectivity is activated at install, MEDRAD Stellant FLEX can automate software updates, including cybersecurity maintenance, and can enable the connected software and services available through the Essential and Comprehensive Tiers.



Everything from the Base Solution:

- + Installation*
- + Initial Clinical Applications Training
- + Remote Connectivity
- + Critical Cybersecurity Patching
- + 12-Month Manufacturer's Warranty

*Installation fees will be listed separately from the MEDRAD Stellant FLEX purchase on the quote and invoice.

BASE SOLUTION TIER

ESSENTIAL SERVICE SOLUTION TIER


ESSENTIAL SOFTWARE SOLUTION TIER

COMPREHENSIVE SOLUTION TIER

EXPLORE ALL TIER OPTIONS

Enhance your Suite with Connected Services or Software

The **Essential Solution Plans** are designed to deliver the fundamentals of connected services or software so that you can effectively manage your daily operations. There are two Essential Solution Plans to choose from, depending on your needs.




BASE SOLUTION TIER

ESSENTIAL SERVICE SOLUTION TIER

ESSENTIAL SOFTWARE SOLUTION TIER

COMPREHENSIVE SOLUTION TIER

The **Essential Service Solution Plan** delivers support and maintenance through remote connectivity and connected services and is ideal for suites with high uptime and performance requirements.



*Installation fees will be listed separately from the MEDRAD Stellant FLEX purchase on the quote and invoice.


†System coverage includes injector and mounting system.

‡Based on the average uptime for all units contracted during the reporting period. Metrics include all fully-supported Bayer product lines excluding Intego.

§Based on a business day (8 AM – 5 PM; local time), excluding Bayer observed holidays. Metrics include all fully-supported Bayer product lines excluding Intego.

¶Reflects single-case issue resolution within one visit.

The **Essential Software Solution Plan** offers access to connected workflow software through a subscription that includes integration and ongoing technical support.



*Installation fees will be listed separately from the MEDRAD Stellant FLEX purchase on the quote and invoice.

[†]Point of Care Interface is required.

[‡]Outcomes dependent on software packages selected.



Unlock the Full Potential of the Connected Suite with Advanced Support, Data and Insights

The **Comprehensive Solution Plan** can elevate productivity, performance and insights in your department with integrated software, services, education and data. In addition, Bayer's clinical and technical experts serve as a valuable extension of your team.



BASE SOLUTION TIER


ESSENTIAL SERVICE SOLUTION TIER

ESSENTIAL SOFTWARE SOLUTION TIER

COMPREHENSIVE SOLUTION TIER

EXPLORE ALL TIER OPTIONS

The **Comprehensive Solution Plan** offers access to connected software, services and advanced clinical support.



The Comprehensive Solution Plan:

- + Optimizes your connected MEDRAD Stellant FLEX solution by pairing hardware, software and support with advanced clinical services.
- + Delivers insight into your suite operations and workflow through connected data.
- + Provides access to Bayer's clinical experts for ongoing collaboration and tailored recommendations.

*Installation fees will be listed separately from the MEDRAD Stellant FLEX purchase on the quote and invoice.

Learn more about what is included in the Comprehensive Solution Maintenance Package and the Comprehensive Software Subscription Package.



A Comprehensive Solution Maintenance Package:

A Comprehensive Solution Maintenance Package delivers hardware and software support with coverage tailored to your performance requirements and service delivery strategy.

Choose between two service packages:

Full-Service Comprehensive Hardware Maintenance*†

Self-Service Comprehensive Hardware Maintenance‡

[Contact your Bayer team to learn more about our customizable options](#)

*System coverage includes injector and mounting system.

†Includes extended preventative maintenance hours (6am-9pm Mon-Fri).

‡This package is anticipated to be commercially available in the Fall of 2026.

§ Full Commercial Release expected in 2026.

¶ Point of Care Interface is required.

A Comprehensive Software Subscription Package:

Includes software integration, full technical support and Advanced Clinical Support.

Software Package

Entitlements

Automated Documentation

- + Point of Care and PACS Interface
- + RIS Interface
- + Speech Recognition Interface
- + Barcode Reader

Contrast Optimization[§]

- + Point of Care and PACS Interface
- + P3T Personalized Dosing Software Bundle
- + Workflow Solutions //Insights

OR BUILD A STANDARD SOFTWARE CONFIGURATION

Standard Software Configuration

Available Software License Options

Choose up to 3 licenses from the available software license options

- + Point of Care and PACS Interface
- + Barcode Reader[¶]
- + RIS Interface[¶]
- + Speech Recognition Interface[¶]
- + P3T Personalized Dosing Software Bundle
- + Workflow Solutions //Insights^{§,¶}

PLUS ADVANCED CLINICAL SUPPORT

Advanced Clinical Support entitles you to 16 hours of ongoing clinical support to drive adoption of software features and solution optimization.

BASE SOLUTION TIER

ESSENTIAL SERVICE SOLUTION TIER

ESSENTIAL SOFTWARE SOLUTION TIER


COMPREHENSIVE SOLUTION TIER

EXPLORE ALL TIER OPTIONS >>>



Ready to get started?

Let's select the right FLEX Solution Plan for you.



[Learn more about
Stellant FLEX](#)



[Contact your Bayer Rep
about customizable
FLEX options](#)



Bayer HealthCare LLC
100 Bayer Boulevard
P.O. Box 915
Whippany, NJ 07981
U.S.A.
Phone: +1-412-767-2400
+1-800-633-7231
Fax: +1-412-767-4120



Manufacturer
Imaxeon Pty Ltd
Unit 1, 38-46 South Street
Rydalmere NSW 2116
Australia
Phone: +1-61-2-8845 4999
Fax: +1-61-2 8845 4998
Imaxeon Pty Ltd is a Bayer company

Reference: 1. James JJ, Tennant SL, Contrast-enhanced spectral mammography (CESM), Clinical Radiology (2018), <https://doi.org/10.1016/j.crad.2018.05.005>. 2. Bayer data on file. Date range 10/1/2024-9/30/2025.

Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

The individuals depicted in this material are actors and not actual health care providers or patients.

Bayer, the Bayer Cross, MEDRAD, MEDRAD Stellant FLEX, Stellant FLEX and Cortenic are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

© 2026 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.

More information on radiologysolutions.bayer.com

EXPLORE ALL TIER OPTIONS