

# The latest technology. Within reach.

Get to the Next Level with TechCARE<sup>™</sup> NxT



## Keeping You Up and Running. And Up to Date.

Bayer in Radiology's TechCARE<sup>™</sup> non-obsolesence programs combine the latest hardware and software enhancements with world-class services to help you maximize your radiology suite potential and focus on your patients.

TechCARE NxT combines foundational services entitlements along with optional enhancements to offer value over the lifetime of your products. For added peace of mind, our team of Clinical Performance Specialists delivers on-site and virtual training to help you and your staff leverage the full potential of your suite of products and services.

#### Foundational entitlements:

- > 3-year equipment service agreement program
- > VirtualCARE<sup>®</sup> Remote Support
- > Injector-level software support and maintenance coverage
- > Access to our Clinical Performance Center
- > One additional virtual clinical training session led by our clinical performance specialists
- > Enhanced access to the myRadiologySolutions customer portal





### With TechCARE NxT, you can elevate your imaging workflow and standardize your equipment without waiting for capital budget approval.

By combining upgrades and uptime, TechCARE NxT provides you with best-in-class monitoring and maintenance, and lets you customize your enhancements to meet your department's specific needs and goals.

### Choose two out of the six available enhancements:

**Upgrade to the latest control room workstation, our informatics-ready platform** The Certegra<sup>®</sup> Workstation 3.0 is designed to increase the functionality and control of your MEDRAD<sup>®</sup> Stellant FLEX CT Injection System.



#### Increase accuracy with automated documentation

Our Certegra<sup>®</sup> Point of Care software enables connectivity to your hospital's information system with patient information from the Modality Worklist, to help match the right patient with the appropriate study. And with the outbound to PACS, contrast details are captured and sent automatically, designed to reduce errors associated with manual data entry.



#### Optimize imaging performance with personalized protocols

Weight-based dosing software automates the calculation of individual protocols based on patient characteristics and contrast concentration to personalize CT studies and increased consistency of individualized protocols.

-(୩)-
<u>\</u>

#### Workflow Solutions //Insights - Actionable intelligence for enterprise management

Our next-generation, cloud-based software provides you with clinical and service dashboards for your connected injectors to enable enhanced management of your department, facility, or enterprise goals.



#### Accurate and Automated data with RIS Interface

The interface allows for contrast injection information to be stored in the RIS and supports accurate coding and invoicing of contrast enhanced procedures for billing purposes.



#### Streamline workflow using Speech Recognition Interface

The customizable interface auto populates patient-specific contrast injection details directly into speech recognition systems and the radiologist's report.

## Workflow Solutions //Insights by Bayer

Workflow Solutions //Insights will provide you with Clinical and Service Key Performance Indicator dashboards for connected injectors that enable immediate, actionable decisions to drive benchmarkbased imaging service improvement.\*

Unlock your own enterprise data to visualize:

- > Injection procedure performance
- > Contrast and syringe utilization
- > Asset management data
- > Equipment service details

With your VirtualCARE<sup>®</sup> connection and our Certegra<sup>®</sup> @ Point of Care software, your injection data is sent securely through VirtualCARE, directly to our Workflow Solutions //Insights cloud-based software, powered by the industryleading ThingWorx platform from PTC., Inc. You can acquire VirtualCARE and Point of Care directly, or through your TechCARE NxT agreement.

\*Will require Point of Care and VirtualCARE Remote Support. A 3-year commitment is required and can be canceled within termination policy. Only available for the MEDRAD® Stellant FLEX CT Injection System.

Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

All patient data that appear in this document are fictitious. No actual patient information is shown. The individuals depicted in this presentation are actors and not actual health care providers or patients.

Bayer, the Bayer Cross, MEDRAD, Certegra, MEDRAD Stellant FLEX, Stellant FLEX, TechCARE, and VirtualCARE are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

© 2021-2024 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.



Bayer HealthCare LLC 100 Bayer Boulevard P.O. Box 915 Whippany, NJ 07981 U.S.A Phone: +1-412-767-2400 +1-800-633-7231 Fax: +1-412-767-4120 More information on radiologysolutions.bayer.com

Manufacturer Bayer Medical Care Inc. 1 Bayer Drive Indianola, PA 15051-0780 U.S.A. Phone: +1-412-767-2400 +1-800-633-7231 Fax: +1-412-767-4120 Indicator Perfo