



Your Intelligent, Total Fluid Delivery System

Syringe and Disposable Kits

Designed for Optimal Workflow

- **Quick, Alignment-Free Attachment and Detachment**
Saves time and maintains aseptic conditions
- **Automatic Plunger Sensing System**
Detects and reacts to syringe state
- **Footless Plunger**
Ensures correct piston-plunger alignment every time

Medrad Disposable Features from Bayer

- Ultra-clear barrel permits visualization of air bubbles
- Patented Medrad® FluiDot® Indicators allow operators to detect the presence of fluid
- Prime tube captures fluid overflow to minimize drip & mess during priming

Additional Options

- Sterile spike minimizes handling of contrast bottle
- Transfer sets economically enhance fluid loading efficiency



medrad® Stellant®
CT Injection System

Sx Single Syringe Kits

Filling Method: Spike

Kits



SSS-SPK
Syringe + fill spike
50 per box

+



SSS-LP-60
60" low pressure tubing +
prime tube
50 per box

+



SVTS
Swabbable-valve
transfer set
50 per box

Suggested Add On's



SSS-CTP-SPK
Syringe + 60" low
pressure tubing with
prime tube + fill spike
50 per box

+



SVTS
Swabbable-valve
transfer set
50 per box

Filling Method: Quick Fill Tube

Kits



SSS-QFT
Syringe + quick-fill
tube
50 per box

+



SSS-LP-60
60" low pressure
tubing with prime tube
50 per box

+



SVTS
Swabbable-valve
transfer set
50 per box

Suggested Add On's



SSS-CTP-QFT
Syringe + 60" low-pressure
tubing with prime tube +
quick-fill tube
50 per box

+



SVTS
Swabbable-valve
transfer set
50 per box

D Dual Syringe Kits

Filling Method: Spike

Kits



SDS-CTP-SPK
Dual syringe kit +60" low
pressure T-tubing with
prime tube + 2 fill spikes
20 per box

+



SVTS
Swabbable-valve
transfer set
50 per box

Suggested Add On's

Filling Method: Quick Fill Tube

Kits



SDS-CTP-QFT
Dual syringe kit with 60"
low-pressure T-tubing with
prime tube + quick-fill tube
20 per box

+



SVTS
Swabbable-valve
transfer set
50 per box

Suggested Add On's

How to Order Quality Medrad® Stellant® Disposable Products from Bayer

Our syringe kit selection allows you to choose the optimum configurations for your procedures. Talk to your Bayer Representative for more information.

New! Order disposables online

The myORDERS website makes it easy to order Bayer disposables and keep track of your orders online. We take care of your orders so you can take care of your patients.

myORDERS
Bayer.com



Bayer

Pharmaceuticals Division

Bayer HealthCare LLC
100 Bayer Boulevard
P.O. Box 915
Whippany, NJ 07981
www.radiologysolutions.bayer.com



Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051 USA
Customer Service/Orders
1-800-633-7231
Customer Service Fax
1-412-767-4120

Bayer reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

Bayer, the Bayer Cross, Medrad®, Stellant®, and FluiDot® are trademarks of the Bayer group of companies.
© 2011, 2014, 2016 Bayer

PP-STL-D-US-0028 November 2016